


Automatic portable stadiometer

BSM 370

Convenient BMI measuring device



Automatic portable stadiometer BSM370 provides precision and movement

Japan Weight Scale Patent Association's approval for accurate testing.



► LCD

- Precise height measurement is given by an ergonomically designed touch bar and detailed measurement sensor.
- Wheels are for short distance mobility.
- Sliding-down & foldable system is for convenient relocation.



Brown



Green



Black

Sliding down & foldable mobility design

Convenient and user friendly design with sliding-down & foldable system helping easy installation.

- Compact light weight design, when fully folded: 1.35m, 14kg.
- Installation does not require any tools.



Accessory equipment for varied use of BSM370

It can be used with various equipment for various purposes.



Wheels for portability.



Remote control unit (optional)

BSM370	2014 / 03 / 15 / 10:27:34
• Height	156,5 cm
• Weight	54,5 kg
• BMI	22,3 kg/m ²
InBody	

Thermal printer (optional)

The BSM370 is designed for portability, both indoors and for shipping to outdoors (such as hospitals and schools).

Remote control unit allows examinee convenient procedure without having to move around from measuring to result check up.

Height, weight, BMI results can be checked and printed.

BSM 370 Specifications

Key Specifications

Measurement Item	Height, Weight, BMI
Height Range	90 ~ 200cm (2ft. 11.4in. ~ 6ft. 6.7in.)
Error Range	± 1mm
Weight Range	2 ~ 250kg (4.5~553.3 lbs)


Feature Specifications

Display Type	Customized LCD
Zero Setting	Automatic zero setting during warming up
Measurement Mode	Weight and Height / Weight only / Height only
Height Measurement	A(Measure height after weight measurement) NA(Press ENTER button after measuring weight to measure height.)
Weight Offset	Unit: 0.1kg, from 0.0kg to 5.0kg Unit: 0.1lbs, from 0.0lbs to 9.9lbs
Unit	Available to choose 'kg/cm' and 'lbs/in.'
Portability	Foldable height measuring part Attaching wheels for carrying easily
Optional Device	Remote control device (2m Cable) Thermal Printer

Other Specifications

Adapter	Power Input	AC 100~240V, 50/60Hz, 1.2A
	Power Output	DC 12V, 3.33A
External Interface	RS-232C 1EA	
Printer	Thermal Printer(Optional, the printers recommended by InBody)	
Dimension	When Touch Bar is down: 380(W) × 510(L) × 1350(H): mm 15(W) × 20(L) × 53.1(H): inch	
	When Touch Bar is up : 380(W) × 510(L) × 2126(H): mm 15(W) × 20(L) × 83.7(H): inch	
Machine Weight	15.5kg (34lbs)	
Measurement Duration	Weight, Height Measurement: 7 sec., Weight Measurement: 2 sec., Height Measurement: 5 sec. (standby time 2 sec.)	
Operation Environment	10 ~ 40□(50 ~ 104□), 30 ~ 75%RH, 70 ~ 106kPa	
Storage Environment	-20 ~ 70□(-4 ~ 158□), 10 ~ 95%RH, 50 ~ 106kPa(No condensation)	

* Specifications are subject to be changed without prior notice.



InBody is a manufacturer of medical measuring and analyzing devices that has acquired over 80 patent rights across the globe.



CE 0120



U.S. patent U.S. 5720296



Canada patent C.N. 2225184



Japan patent



ISO13485



ISO9001



Korea Food & Drug Administration

InBody

InBody Co., Ltd. [HEAD OFFICE]

TEL: +82-2-501-3939

FAX: +82-2-578-2716

Website: <http://www.inbody.com>

E-mail: info@inbody.com

InBody [USA]

TEL: +1-323-932-6503

FAX: +1-323-952-5009

Website: <http://www.inbodyusa.com>

E-mail: info@inbodyusa.com

InBody Japan Inc. [JAPAN]

TEL: +81-03-5298-7667

FAX: +81-03-5298-7668

Website: <http://www.inbody.co.jp>

E-mail: inbody@inbody.co.jp

Biospace China. [CHINA]

TEL: +86-21-64439738, 9739, 9705

FAX: +86-21-64439706

Website: <http://www.biospacechina.com>

E-mail: info@biospacechina.com